

## Metabolomics Extraction from Saliva Samples

1. Time of sample collection must be justified to decrease the interference from the food and drink. It is better to fasten the animals for 1-2hr.
2. Visual inspection must reveal no contaminations, e.g. hair and blood, in the samples.
3. The samples must be stored in a freezer at  $-80^{\circ}\text{C}$  until extraction.
4. Equal samples amount must be used during the extraction.
5. Cold extraction buffer with thawing and freezing technique

- Prepare extraction solvent: Methanol: dichloromethane: ethyl acetate 1:2:3.
- Extraction solvent kept in  $-80^{\circ}\text{C}$  freezer for at least 2hr.
- 200  $\mu\text{l}$  of the sample is mixed with 900  $\mu\text{l}$  of pre-cooled extraction solvent in an Eppendorf tube.
- Vigorous mixing is done using vortex for 1 minutes.
- The tube is transferred into dry ice or freezer ( $-80^{\circ}\text{C}$ ) for 30 min.
- The sample is thawed in an ice bath then it is sonicated for 5 min at 20-30 kHz. (temperature must be not greater than  $20^{\circ}\text{C}$ )
- The sample is centrifuged at 10,000 for 10 min  $4^{\circ}\text{C}$ .
- The supernatant is transferred to a new tube.
- The sample is centrifuged at 10,000 for 10 min. The supernatant is transferred to a new tube.
- The supernatant evaporated using vacuum rotary evaporator at  $30^{\circ}\text{C}$  to approximately 25  $\mu\text{l}$  volume.
- The sample reconstituted in 200 $\mu\text{l}$  solvent (water: methanol: acetonitrile 2:1:1).
- Centrifuged at 10,000 for 5min. The supernatant is transferred to analysis tube.
- Analysis, done using 5-25  $\mu\text{l}$  injection volume (depending on the complicity of the sample), is applied to the LC-MS/MS. This will be used for positive and negative mode.

Proteomics & Metabolomics Unit  
[proteomics.lab@57357.org](mailto:proteomics.lab@57357.org)

Children's Cancer Hospital 57357  
Cairo, Egypt

<https://www.57357.org/proteomics-unit/>

- Mobile phase depends on the mode:

Positive mode mobile phase	Negative mode phase
DI-Water contains 0.1% FA	5mM ammonium format containing 1% methanol. Then the pH will be adjusted to 8 using ammonium hydroxide.

Proteomics & Metabolomics Unit  
[proteomics.lab@57357.org](mailto:proteomics.lab@57357.org)

Children's Cancer Hospital 57357  
 Cairo, Egypt

<https://www.57357.org/proteomics-unit/>

- Quality control (QC) samples, which are a mixture of equal volume taken from each real sample, also underwent LC-MS/MS analysis for quality assurance of the experiment.
- It is better to spike the sample with internal standards. These internal standard is better that not included in the endogenous metabolites. You can also use another one depending on experimental design.