Metabolomics Extraction from Stool Samples

- 1. Store collected fecal samples at -80°C until usage.
- 2. Thaw samples on ice.
- 3. Take 100 mg from each sample in new Eppendorf.
- 4. First Extraction:
 - Add 400 µl of extraction solvent (acetonitrile: methanol, 3:1) to the samples.
 - Vortex for 2 minutes.
 - Sonicate for 5 minutes.
 - Centrifuge for 10 minutes at 4°C at 10,000 rpm.
 - Transfer the supernatant to new Eppendorf.
- 5. Second Extraction: (Repeat the above steps)
 - Add 400 μ l of the extraction solvent to the pellet and homogenize.
 - Sonicate for 15 minutes.
 - Centrifuge for 10 minutes at 4°C at 10,000 rpm.
 - Transfer the supernatant to the Eppendorf of the first extraction step.
- 6. Speed-vac the supernatant at 30° C.
- 7. Reconstitute in reconstitution buffer (water: methanol: acetonitrile, 2:1:1) at ratio 1mg/1ml.



Proteomics & Metabolomics Unit proteomics.lab@57357.org

Children's Cancer Hospital 57357 Cairo, Egypt

https://www.57357.org/proteo mics-unit/