

Stable isotope dimethyl labeling

[The protocol is designed to label 25-30 µg tryptic peptides]

[All preparation should be fresh]

1. Extract proteins from cell lines or biological tissue using 8M urea in 100 mM TEAB, protease inhibitor cocktail.
2. Quantify proteins using BCA method.
3. Start with **25- 30 µg** protein extract, proceed with amine- free, In- solution trypsin digestion protocol as follows;

#scale up or down based on volume (optimally 30-50 µl)

- a) **30 µl** 8 M Urea buffer [240 mg urea, 100 µl of 100 mM TEAB, 220 µl MilliQ, protease inhibitors]
- b) **0.2µl** 1M TCEP (tris (2-carboxyethyl) phosphine) in milliQ
mol. Mass 250.19 g mol⁻¹
- c) Incubate at room temperature for 30 minutes
- d) **0.5 µl** of 500 mM IAA in milliQ
- e) Incubate at room temperature for 30 min in dark
- f) **90 µl** 100 mM TEAB to get final urea dilution around 2M
- g) Trypsin (1:20): prepare 0.5ug/µl trypsin solution and add 2 µl
- h) Add **4µl** 100 mM CaCl₂
- i) Incubate overnight at 37°C in shaker (600 rpm)
- j) Spin down for 30 min at max speed at room temperature

Don't acidify the sample

4. Start with trypsinized peptides in TEAB and urea buffer without final acidification.

[above; total 125 µl contains 30 µg with a final conc 0.24 ug/µl]

Proteomics & Metabolomics Unit
proteomics.lab@57357.org

Children's Cancer Hospital 57357
Cairo, Egypt

<https://www.57357.org/proteomics-unit/>

5. Take **100 µl** [contains 24 µg digested peptide] for subsequent labeling, make sure pH 5.0~8.0. (check with pH paper)
6. Add **4µl** of 4% (vol/vol) of either formaldehyde (CH₂O), Deuterated formaldehyde (CD₂O), or C13, deuterated formaldehyde (¹³CD₂O) to be labeled with light, intermediate, or heavy dimethyl, respectively. (CH₂O is diluted **1:8** (2µl F+ 16 µl milliQ), CD₂O and ¹³CD₂O are diluted 1:4 (2µl F+ 8 µl milliQ))
7. Mix briefly, spin down the sample
8. Prepare cyanoborohydride, cyanoborodeuterate solutions to 0.6M (NaBH₃CN, and NaBD₃CN) (prepare fresh). [0.018 g/500 µl milliQ]
9. Add **4 µl** of 0.6M NaBH₃CN to light and intermediate labeling, and **4 ml** of 0.6 M NaBD₃CN heavy labels

Proteomics & Metabolomics Unit
proteomics.lab@57357.org

Children's Cancer Hospital 57357
 Cairo, Egypt

<https://www.57357.org/proteomics-unit/>

Light	Medium	Heavy
+28.0313	+32.0564	36.0757
CH ₂ O	CD ₂ O	¹³ CD ₂ O
NaBH ₃ CN	NaBH ₃ CN	NaBD ₃ CN

10. Incubate at RT for 1 hour, (preferably in the hood in a shaker)
11. Quench reaction (in the hood) by the addition of **15µl** of 200mM ammonium bicarbonate (or 1% ammonia solution)-mix [7.9 g/500 µl]
12. Add 8µl of formic acid to further quench reaction and acidify sample.
13. Stage tip the sample (30 ug/130 µl with conc 0.2 ug/µl)